

Capsule	Typical Fill Weights (gm) Powder Density <sup>a</sup>			Body Volume⁵	Locked Length		External Diameter <sup>c</sup>	Cut Length	Single Wall Thickness	Dome Thickness	Weight (Avg. of 100)
	0.45 (light)	0.7 (typical)	1.00 (heavy)	(ml)	± 1.0 (mm)		± 0.10 (mm)	± 0.6 (mm)	± 0.02 (mm)	± 0.05 (mm)	± 10% (gm)
602340	12.6	19.6	28.0	28.0	87.6	Сар	23.33	30.0	0.22	0.29	2.22
10pk	.2.0	10.0	20.0	20.0	07.10	Body	22.72	76.8	0.22	0.30	
602343 24pk	3.4	5.3	7.5	7.5	57.0	Cap	15.27	21.0	0.20	0.24	0.85
						Body	14.67	51.0	0.19	0.24	

<sup>&</sup>lt;sup>a</sup> Actual Fill Weights vary and depend on powder characteristics and tamping/compression during filling.
<sup>b</sup> Capsule volume in ml is of the capsule body only and does not include the cap dome.
<sup>c</sup> Capsule body and cap have a slight taper and diameter given is the maximum diameter at the open end.